



## 10W 12 Needles Microneedle Therapy Machine For Skin Rejuvenation

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: WZ-3
- Minimum Order Quantity: 1
- Price: USD/60-100
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



#### Product Specification

- Power Consumption: 10W
- Color: White
- Function: Skin Rejuvenation, Scar Removal, Wrinkle Removal
- Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual
- Needle Length: 0.25mm-2.5mm
- Needle Number: 12
- Power Supply: AC 100-240V, 50/60Hz
- Usage: Home, Salon
- Name: Microneedle Therapy Machine
- Highlight: **10W Microneedle Therapy Machine, 12 Needles Microneedle Therapy Machine, Skin Rejuvenation Microneedle Therapy**



#### More Images





## Product Description

### 10W White Microneedle Machine For Skin Rejuvenation With 12 Needles

#### Product Description:

The Microneedle Therapy Machine features 5 speed levels, allowing you to customize the intensity of your treatment. Whether you prefer a gentle massage or a more intense experience, the machine has you covered. The power consumption of the machine is only 10W, making it energy-efficient and environmentally friendly.

The Microneedle Therapy Machine is an excellent choice for those looking to improve the appearance and texture of their skin. The 12 needles on the machine penetrate the skin, stimulating collagen production and promoting cell regeneration. This can result in firmer, smoother, and more youthful-looking skin.

In addition to its skin rejuvenation benefits, the Microneedle Therapy Machine is also a great tool for relaxation. The electric wand massager provides a soothing massage that can relieve tension and reduce stress. Use it on your face or body to unwind after a long day or to help you sleep more soundly at night.

The Microneedle Therapy Machine is easy to use and comes with clear instructions. Simply apply your favorite serum or moisturizer to your face or body, adjust the speed level to your desired intensity, and move the machine over the treatment area. The machine is lightweight and ergonomically designed, making it comfortable to hold and use for extended periods of time.

Overall, the Microneedle Therapy Machine is an excellent choice for anyone looking for an effective and versatile electric wand massager. With its 12 needle design, 5 speed levels, and ability to treat both the face and body, it is a must-have for those looking to improve the appearance and texture of their skin while also enjoying the benefits of relaxation and stress relief.

#### Features:

Product Name: Microneedle Therapy Machine

Needle Length: 0.25mm-2.5mm

Needle Number: 12

Usage: Home, Salon

Power Supply: AC 100-240V, 50/60Hz

Model: MTM-01

BIO Microcurrent Machine

Hot Cold Therapy Machine

Butt Vacuum Therapy Machine

#### Technical Parameters:

Needle Number	12
Package Contents	Microneedle Therapy Machine, Power Adapter, User Manual
Usage	Home, Salon
Power Consumption	10W
Power Supply	AC 100-240V, 50/60Hz
Material	Stainless Steel
Speed Levels	5
Color	White
Treatment Area	Face, Body
Needle Length	0.25mm-2.5mm

#### Applications:

One of the primary uses of the Microneedle Therapy Machine is as an electric wand massager. With its five-speed levels, this machine can be customized to provide the perfect amount of pressure and stimulation for any user. It is great for relieving muscle tension and promoting relaxation.

Another application for the Microneedle Therapy Machine is as a BIO microcurrent machine. This technology uses low-level electrical currents to stimulate the skin and promote collagen production. This helps to reduce the appearance of fine lines and wrinkles, as well as improving skin texture and tone.

In addition to its use as an electric wand massager and BIO microcurrent machine, the Microneedle Therapy Machine can also be used

for microneedling. This technique involves using tiny needles to create micro-injuries in the skin, which then stimulates the body's healing response and promotes collagen production. This can help to improve the appearance of scars, stretch marks, and other skin imperfections.

The Microneedle Therapy Machine is suitable for use in a variety of settings, including spas, salons, and at home. It is powered by an AC 100-240V, 50/60Hz power supply, making it easy to use anywhere in the world.

Overall, the Microneedle Therapy Machine is a versatile and effective device that can be used for a variety of purposes. Whether you are looking to relax your muscles, improve your skin, or reduce the appearance of scars, this machine is a great choice.

### Customization:

Model Number: WZ-3

Place of Origin: China Guangzhou

Power Consumption: 10W

Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual

Usage: Home, Salon


Color: White

Model: MTM-01

Enhance your Microneedle Therapy Machine with our BIO Microcurrent Machine, gold RF microneedle, and Butt Vacuum Therapy Machine options. Make it stand out and suit your needs perfectly!



## Product Parameters



### ● Specification

Brand	Hydra pen	Model	H5
EMS intensity level	4 levels	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② Inductive base ③ Hydra pen H5 ④ Cable ⑤ Adapter ⑥ Manual ⑦ Retaining ring* 2 pcs ⑧ 12 pins * 2 pcs

### Adjustable liquid output

please adjust the discharge of liquid base on the consistency of liquid.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebain Village, Baiyun District, Guangzhou