



Home Unwanted Hair Removal Laser Machine 1200W With Pulse Duration 10-400 Ms

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: TM-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Unwanted Hair Removal Laser Machine
- Pulse Duration: 10-400 Ms
- Machine Power: 1200W
- Energy Density: Up To 120 J/cm²
- Display: 15-inch Touch Screen
- Machine App: Permanent Laser Hair Removal Machines
- Style: Stationary, Vertical
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Coolingtype: TEC+Sapphire Contact Cooling
- Highlight: Home Unwanted Hair Removal Laser Machine, Unwanted Hair Removal Laser Machine 1200W, CE Unwanted Hair Removal Laser Machine



More Images



Product Description

Professional Unwanted Hair Remover Device With Pulse Duration 10-400 Ms

Product Description:

Equipped with a powerful 600W laser, this Laser Hair Treatment Device offers efficient and effective hair removal for all skin types. The TEC+Sapphire Contact Cooling system ensures a comfortable and painless experience by cooling the skin during the treatment process. Additionally, the machine is equipped with 2pcs Fans for added cooling and safety.

Our Laser Hair Elimination Machine is designed with safety in mind. The Skin Contact Cooling feature ensures that the skin is protected from the heat generated by the laser, while the Automatic Power Cut-off feature ensures that the machine automatically shuts off when it detects any abnormality in the treatment area.

With a Machine Power of 1200W, this Laser Hair Treatment Device provides fast and efficient hair removal, making it ideal for busy individuals who want to achieve smooth, hair-free skin without having to spend hours in a salon or spa. The device is easy to use and comes with detailed instructions for hassle-free operation.

Experience the convenience and effectiveness of our Laser Hair Elimination Machine, the perfect solution for anyone looking to remove unwanted hair in the comfort of their own home. Say goodbye to traditional hair removal methods and hello to smooth, hair-free skin with this innovative and highly effective Laser Hair Treatment Device.

Features:

Product Name: Hair Removal Laser Machine

Style: Stationary, Vertical

Energy Density: Up To 120 J/cm²

Cooling Fan: 2pcs Fans

Machine App: Permanent Laser Hair Removal Machines

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

Introducing the Advanced Hair Removal System - the ultimate Hair Depilator Machine for achieving smooth skin. This powerful laser machine boasts a stationary, vertical design and an energy density of up to 120 J/cm². With 2 cooling fans and safety features like skin contact cooling and automatic power cut-off, this Smooth Skin Laser Machine is both safe and effective. Use the machine app for easy, permanent laser hair removal from the comfort of your own home.

INTRODUCTION OF ACCESSORIES





INTERFACE

//



PRODUCT DETAILS

Learn more about the product



APPLICATIONS



**Hair
Removal**




**Hemangioma
Removal**



**Blood Vein
Removal**



**Vascular
Removal**



**Skin
Rejuvenation**



**Cure Ringworm
of The Nails**

Technical Parameters:

Machine Power	1200W
Machine App	Permanent Laser Hair Removal Machines
Pulse Duration	10-400 Ms
Cooling type	TEC+Sapphire Contact Cooling
Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Power	600W
Skin Types	I-VI
Style	Stationary, Vertical
Cooling Fan	2pcs Fans
Energy Density	Up To 120 J/cm2

This Laser Hair Treatment Device is a Painless Hair Removal System that utilizes a Permanent Hair Removal Device. It is suitable for skin types I-VI and features a cooling system with TEC+Sapphire Contact Cooling and 2pcs Fans. It also includes safety features such as Skin Contact Cooling and Automatic Power Cut-off.

Applications:

The LZR TM-1 hair removal laser machine is equipped with a 15-inch touch screen display, making it easy to operate and control. Users can easily navigate through the settings and select the desired treatment mode. The machine also has two cooling fans that help to keep the device at a stable temperature during use, preventing overheating and ensuring maximum efficiency. The cooling type used is TEC+Sapphire Contact Cooling, which provides a comfortable and painless experience for the user.

This professional hair removal device is perfect for use in various occasions and scenarios. It is suitable for use in beauty salons, spas, and clinics, where hair removal treatments are offered. The LZR TM-1 can be used to remove hair from various parts of the body, including the legs, arms, underarms, bikini area, and face. It is a reliable and efficient unwanted hair remover device that can be used on all skin types, including sensitive skin.

The LZR TM-1 hair removal laser machine is also perfect for personal use. It can be used at home by individuals who want to remove

unwanted hair from their body. The device is easy to operate, and users can adjust the settings to their desired level. This hair removal device saves time and money, as users no longer need to visit a beauty salon or clinic for hair removal treatments. In conclusion, the LZR TM-1 hair removal laser machine is a professional and efficient hair depilator machine that can be used in various occasions and scenarios. It is perfect for use in beauty salons, spas, and clinics, as well as for personal use at home. The device uses the latest laser technology, has a touch screen display, and is equipped with two cooling fans that ensure maximum efficiency and comfort. Get your LZR TM-1 hair removal laser machine today and say goodbye to unwanted hair for good!

Customization:

Support and Services:

The Hair Removal Laser Machine is a high-quality product that requires minimal maintenance. However, in the rare event that you experience technical difficulties with your machine, our technical support team is available to assist you.

Our technical support services include:

Phone and email support

Remote diagnostic testing

On-site technical support (additional charges may apply)

In addition to our technical support services, we also offer product maintenance and repair services. These services include:

Regular maintenance checks to ensure optimal product performance

Product repair services for any damages

Product upgrades and modifications (additional charges may apply)

Our ultimate goal is to ensure that you are completely satisfied with your Hair Removal Laser Machine product and its performance. If you have any questions or concerns, please do not hesitate to contact our technical support team.

FAQ:

Q: What is the brand name of this hair removal laser machine?

A: The brand name of this hair removal laser machine is LZR.

Q: What is the model number of this hair removal laser machine?

A: The model number of this hair removal laser machine is TM-1.

Q: Where is this hair removal laser machine made?

A: This hair removal laser machine is made in China Guangzhou.

Q: Can this hair removal laser machine be used on all skin types?

A: This hair removal laser machine is suitable for all skin types except for very dark skin.

Q: How long does it take to see results from using this hair removal laser machine?

A: It typically takes 4-6 weeks to see results from using this hair removal laser machine.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaity-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou