



Continuous / Fractional Ultra Pulse CO2 Laser Machine For Skin Resurfacing

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: Q9
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: RF Fractional CO2 Laser
- Laser Arm Origin: From South Korea
- Power: 40 Watts
- Operation System: Ultrasonic
- Laser Working Mode: Continuous
- Aiming Beam: 650nm, ≤5mw
- Output Mode: Ultra Pulse, Single Pulse, Fractional
- Laser Type: Fractional CO2 Laser
- Language: Multi-language
- Highlight: Continuous Ultra Pulse CO2 Laser, Fractional Ultra Pulse CO2 Laser, Ultra Pulse CO2 Laser Machine



More Images



Product Description

40 Watt Continuous / Fractional Ultra Pulse CO2 Laser Machine For Skin Resurfacing

Product Description:

This machine operates on an Ultrasonic system, which makes it incredibly precise and effective. It also features three different output modes: Ultra Pulse, Single Pulse, and Fractional. This means that you can choose the mode that best suits your needs and preferences. One of the key features of this machine is its Fractional Working Mode. This mode allows for interval continuous working, which means that you can treat larger areas of skin without having to stop and start the machine repeatedly. This makes the treatment process much more efficient and convenient.

As a Fractional CO2 Laser, this machine is capable of delivering incredibly precise and powerful laser energy to the skin. This energy helps to stimulate the skin's natural healing processes, which can lead to smoother, firmer, and more youthful-looking skin. In addition to its Fractional CO2 Laser capabilities, this machine also features Elight Rf Laser technology. This technology helps to further enhance the results of the treatment by delivering targeted energy to the deeper layers of the skin. Overall, the Skin Rejuvenation Fractional CO2 Laser Machine is an excellent choice for anyone looking to improve the appearance of their skin. With its powerful output, precise operation, and advanced technologies, this machine is sure to deliver impressive results.

Features:

Product Name: RF Fractional CO2 Laser

Fractional Working Mode: Interval Continuous Working

Net Weight: 45kg

Power: 40 Watts

Style: Portable

Language: Multi-language

Other keywords: Bystronic CO2 Laser, Medical Fractional CO2 Laser Machine, IPL RF ND YAG Laser



Technical Parameters:

Laser Arm Origin	South Korea
Laser Type	Fractional CO2 Laser

Fractional Working Mode	Interval Continuous Working
Size	360x450x1250mm
Light Guiding System	7 Joints Arm
Power	40 Watts
Laser Working Mode	Continuous
Style	Portable
Language	Multi-language
Headpiece	Fractional, Vaginal, Normal Cut

Applications:

The Elight Rf Laser Machine is perfect for various skin concerns. It can be used for skin resurfacing, wrinkle reduction, scar removal, and stretch mark reduction. The machine is also effective in treating sun damage, age spots, and acne scars. It uses a high-powered CO2 laser that penetrates deep into the skin, stimulating collagen production and promoting skin rejuvenation.

The CO2 Fractional Machine is ideal for individuals who want to improve their skin's texture and tone. It can be used to treat a variety of skin concerns, including hyperpigmentation, wrinkles, and acne scars. The machine works by creating tiny holes in the skin, which triggers the body's natural healing process. This process promotes collagen production, which helps to improve skin texture, reduce wrinkles, and minimize pores.

The RF Fractional CO2 Laser is perfect for beauty clinics and spas that want to offer a wide range of skin treatments. It is also ideal for individuals who want to improve their skin without undergoing invasive procedures. This machine is easy to use, and its multi-language options make it accessible to people from different regions. Whether you are looking to reduce wrinkles, improve skin texture, or treat acne scars, the Elight Rf Laser Machine is an excellent choice.

Customization:

FAQ:

Q: What is the place of origin of the RF Fractional CO2 Laser?

A: The RF Fractional CO2 Laser is made in China Guangzhou.

Q: What is the purpose of the RF Fractional CO2 Laser?

A: The RF Fractional CO2 Laser is used for skin resurfacing, scar removal, wrinkle removal, and skin rejuvenation.

Q: How does the RF Fractional CO2 Laser work?

A: The RF Fractional CO2 Laser uses laser technology to create tiny holes in the skin. This triggers the body's natural healing process, which promotes the growth of new, healthy skin cells. Additionally, the RF technology stimulates collagen production, which helps to improve skin texture and reduce the appearance of wrinkles.

Q: Is the RF Fractional CO2 Laser painful?

A: The RF Fractional CO2 Laser may cause some discomfort during the treatment, but most patients report only mild pain. A topical anesthetic can be applied prior to the treatment to help minimize discomfort.

Q: How many treatments are required to see results?

A: The number of treatments required depends on the individual patient and their specific skin concerns. Generally, patients will see results after just one treatment, but multiple treatments may be necessary for optimal results.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou